



Dry vacuum systems

Tru-Vac single TV-12

Specifications



Capacity: Up to 6 Users

Power: 2.4 HP

Warranty: 10 Years

Voltage: 115 / 230 V AC*

Flow Rate: 85 CFM

Current: 20 / 10 Amps

Noise Level: 61 dB

Tank size: 40" Tall / 12 Gal.

Vacuum ("Hg): 8.5" Hg

Tank Weight: 43 lbs

Chassis Footprint: 21"x 20"x 16"

Completely Waterless

Chassis Weight: 85 lbs

Completely Oil-Less

Please specify voltage requirement when placing order