

Made in the
USA!



The **EAGLE VORTEX** line of Oil Free Dental Compressors

Features

- Oil-less design
- High capacity, permanently lubricated and sealed bearings
- All aluminum design
- Automatic Thermal Protection
- Enclosed crankcase
- Low operating sound level
- Long, service free life
- Increased field service capability
- High volumetric efficiency
- 90 degree V-design-low vibration, high efficiency rating, & full rated pressure restart capability



Direct fit on any EGL- 12
Models, MSV Models,
and even Tech West,
DCI and more!

Heads up!

Vortex Reciprocating Piston type compressor heads like the ones in these images. These very reliable pumps are paired with a Baldor Electric motor and have a very long service life. The revolutionary Vortex Series of articulating piston air compressors raises the bar for innovation and durability.

