

# SPECIFICATIONS



## SITE REQUIREMENTS

Voltage	115 Volts AC
Circuit breaker rating	20 Amp
Minimum wire size	14 Gauge
Minimum trunk line size	½"
Max ambient temp.	100° F
Minimum ambient temp.	50° F

This equipment must be installed in a well ventilated mechanical room to prevent overheating and failure not covered by the warranty.

Description	Specification
Model #	EGL-T6
Max users	4
Air dryer	Yes/Dessicant/Auto purge
Number of heads	2
Total horsepower	1.5
CFM	5.05
Tank gallons	8 g
Voltage	115 V AC
Amps	11
Decibels	67
Weight	110 lbs.
Height	28"
Width	18"
Depth	24"